

# CepaCom

Mobile Communicator



*The founders of Cepa Mobility (2003): Niels, Tarik, Per og Thomas*

## Features:



Speech Machine



Media Player



Mobile Phone



Text Messaging

**Cepa Mobility®**  
- Enabling the Disabled -

[www.cepamobility.com](http://www.cepamobility.com)

# CepaCom


Enables You ...

CepaCom is a small **plug-and-play** mobile communicator. It mounts on to your wheelchair or bed easily and enables you to communicate with your family, helpers and friends without additional help. Choose between...

- **Mobile Phone** and **Text Message** functionalities. Efficient and independent communication. Get involved in work and other social activities.
- **Speech Machine** with high quality text to speech. CepaCom can read your everyday sentences and text messages out loud.
- **Media Player**. Listen to e-books, interviews, music, lectures, meetings and learn more.



#### Speech Machine Example:



14:22 - *I would like to go out for dinner*

CepaCom is easy to operate with a 01-switch, a joystick\* or directly by hand. Due to the wireless navigation you can easily move the CepaCom between a wheelchair and a bed. The CepaCom is small and discreet. It can be charged directly from your wheelchair, a wall plug or mobile power pack.

The cables, plug-ins and the Pocket PC are well protected from dust, rain and strong knocks or bumps when mounted on a wheelchair.

# What do our customers say? ...

## Like having my voice back

“Cepa Mobility’s communicator is incredibly easy to use. It s helped me keep in touch with my network. With the speech-machine I have an easier time communicating with assistants and others in my surroundings. It has been like having my voice back.”

Jens Maahr (42) ALS patient

## Improves my work

“Now I can do my job like other people. As a nutritionist and therapist I have sworn professional secrecy. Therefore it is of great importance that I am able to communicate with my clients over the phone and through text messages without help from others. It enables me to keep my work private and not involve my helpers. With the CepaCom I am also a more independent person with more privacy than ever before. I can now send my helpers away when I want to be alone and I know that with the CepaCom I can just call them when I need them.”

Mette (24) CP

## Kind and friendly service

“When you slowly lose your ability to move and speak due to ALS, CepaCom enables you to stay in contact with your friends and family – it is fantastic. Unfortunately my friend has this terrible disease, but in spite of that he is able to receive, write and send text messages from his hospital bed by using only his toe to operate the CepaCom. Being able to receive and answer the many greetings from the outside world means everything to him. It is also nice for family and friends to be able to keep in contact when he is away from the hospital. Thank you for your fast and kind service.”


Maiken (42) Relative to ALS patient

## We can communicate easier and faster

“CepaCom contributes to a life independant of the many nurses surrounding us. Independance means a lot to me as a spouse, since it affects Jens in a positive way and heightens his dignity as a human being. So we have both been helped a lot in our daily lives with Cepa Mobilitys solution. Im no longer powerless when we communicate and I become less angry. We have Cepa to thank for that. We have both experienced a very fine service incredible drive to make things work.”

Heidi Storbank (42) Wife to ALS patient

## Text Message Example:



16:43 - Writing a message to Sam

## Media Player Example::



16:22 - This is an exciting book

**Pedagogical values:** Better, private and independant communication  
• Better understanding of your own self • Clearer demarcation of your own identity • Maintain and expand your network • Pass on your own stories and experiences • Breaking down barriers between disabled and non-disabled.

CepaCom is very easy to operate and has enlarged icons and text. You can setup and adjust CepaCom to your personal needs and change them as you improve your operating skills. CepaCom has an intelligent dictionary which helps you write text messages. Go to [www.cepamobility.com](http://www.cepamobility.com) for more information.

# CepaCom

The New Mobile Communicator

## The People and the Story

Cepa Mobility was founded in the year 2003 when a young student named Niels Martensen started working as an assistant for Tarik Ghiyati who has cerebral palsy. Tarik had a strong wish to use a mobile device on his own. Together they started developing the CepaCom. They were later joined by Per Rasmussen and Thomas Pedersen. Per is a software developer and Thomas is a hardware engineer. Today they run Cepa Mobility together and by using the latest mobile and wireless technologies they continue to develop communication solutions for people with special needs.

## Customer Focus

We have a continuous dialogue with our customers when we test and develop new features. In this way we get the best results and make sure that our solutions are optimised for the user and the surrounding environment.

## How do I get CepaCom?

CepaCom is an AAC device and can be subsidised as such. Contact your country specific distributor for prices and information about how to get a CepaCom.

## Try out CepaCom

We are happy to demonstrate CepaCom for you. Contact your country specific distributor for information on how to get a demonstration.

## How do I Navigate CepaCom?

**01-Switch:** The different menus roll over the screen in a scanning process. When the menu item you want to select is situated in the middle section of the screen you activate it by pressing the 01-switch. This way you can navigate the different menus in the CepaCom.

**Joystick:** It is also possible to navigate the CepaCom by using the joystick on your wheelchair. Navigation is faster this way compared to switch interaction. The following wheelchair support CepaCom: RolTec Vision+.

## Mounting

CepaCom is mounted in an ergonomically correct manner by professionals. Our distributors have many years of mounting experience and they are familiar with protecting CepaCom from the impacts of daily life.

## Prices and Specifications

Contact us directly or find your country specific distributor for further information at [www.cepamobility.com](http://www.cepamobility.com).

### Distributor

### Main Office

Cepa Mobility Aps™  
Njalsgade 88 3.3 101  
2300 Copenhagen S  
Denmark

Sales: +45 70 27 74 20

[www.cepamobility.com](http://www.cepamobility.com)

Mail: [info@cepa.dk](mailto:info@cepa.dk)